


SYMEDRIX

Aortic Ultra[®]

AAA WIRE TECHNOLOGY



AORTIC ULTRA Guidewire

Symedrix' AORTIC ULTRA is specifically designed for aortic interventions.

Designed for aortic interventions

The AORTIC ULTRA is a coated ultra stiff .035" guide wire, which is developed for aortic interventions. The PTFE coating enables a smooth delivery, while the ultra stiffness of the AORTIC ULTRA provides excellent support to maneuver catheter-based systems securely to the desired location.

Ultra stiff main body


The AORTIC ULTRA's unique ultra stiff main body provides excellent guidance and support even in challenging anatomies for aortic interventions.

Coating

For the AORTIC ULTRA a special kind of coating is used. THE PTFE coating is thicker than usual. Due to the enhanced number of PTFE particles, the AORTIC ULTRA ensures a more frictionless and safe delivery.

Additional Features

- Atraumatic Tip for safe positioning
- Tip visibility fo precise positioning
- Safe kink-free transition between ultra stiff main body and tip



Reference
Testing performed by Staatliche Materialprüfungsanstalt Darmstadt; test data on file at Symedrix GmbH

AORTIC ULTRA 0.035" GUIDE WIRE

REF Number*	Diameter	Length	Tip Curve Diameter
AU-035-260-C	0,035 inch	260 cm	20 mm J-Tip

0,035 inch / 0,89 mm

Outer Diameter

J-Tip

TCD
20 mm

*(Box of 5 units)

SYMEDRIX

SYMEDRIX GmbH

Bahnhofstrasse 32

D-82041 Oberhaching

www.symedrix.com

Tel: +49 89 666 58 79 29

Fax: +49 89 54 85 80 10

info@symedrix.com